



## ALL-IN-ONE DELIVERY

Scan the QR code to learn more about the unique technology behind delivering an all-in-one delivery tray of veneers at the same time.



FIRST FIT

Powered by HENRY SCHEIN

## LOCATE A LAB NEAR YOU

Scan the QR code to view the list of FirstFit Authorized Lab Providers.



## FOR FREE INFORMATION:

800.921.4806  
ext. 131  
www.dps.  
li/a/68J-149

# Single-Visit Veneers

FirstFit technology allows you to prep and seat restorations in the same appointment

It's been said that a cosmetic dentist is part artist and part clinician, especially when it comes to tasks like building direct composite veneers. But advances in technology are transforming the cosmetic dentistry landscape with revolutionary new products that ensure predictable and consistently esthetic restorations—every time.

FirstFit is one of those solutions. The innovative guided prosthetics delivery system provides the clinician with 3D-printed prep guides and the final, ready-to-seat restoration—veneers, bridges, or crowns—at the same time. In the case of veneers, the clinician receives a simple, all-in-one tray of veneers to be delivered at the same time. The delivery tray provides full visibility of the veneers and teeth so that seating, angulation, and margins are all perfectly aligned.

## Efficiency by 3D-Printed Design

Dr. Ross W. Nash, an international lecturer on esthetic dentistry, sees many patients for veneers at his Charlotte, NC, cosmetic practice. "It takes a long time to put 10 veneers in at one time—almost an hour a tooth," he said. "With FirstFit, you can fit up to 10 veneers and deliver them all in at once. They're in a tray, they're perfectly seated, and you can make sure they're all seated before you light-cure. Then you simply remove the tray. It's very efficient."

The FirstFit guided veneers delivery system, which includes 3D-printed guides, a model, ready-to-seat veneers in a delivery tray, and 2 burs, minimizes tooth preparation so the teeth are reshaped to exact specifications for the veneers—reducing finishing time and minimizing patient discomfort.

"For a dentist who is looking to place more all-ceramic veneers, FirstFit is a good tool to promote that in their practice," said Dr. Nash. "It's a very straightforward system, and I think it's something that dentistry needs."



Ross W. Nash,  
DDS

"It's a very straightforward system, and I think it's something that dentistry needs."

## The Patient Experience

While FirstFit offers clinicians like Dr. Nash a simplified and predictable technique for placing veneers, patients love the fact that they can restore their smiles in a single visit. "Patients are ecstatic when then can come in, and in one day, get all-ceramic veneers that usually take 2 appointments several weeks apart," said Dr. Nash. "My patients who I've placed these on have been very happy."

Shaharyar Fiza, DMD, who practices with Dr. Nash, agrees. "I like the single appointment aspect of FirstFit and the convenience of the guided system. And I know patients appreciate not having to wear temporaries," he said.

